**Symbols / Symbolerklärungen / Simbolos y su significado /
Explications des symboles / Spiegazione dei simboli**

Conformity mark
Konformitätskennzeichen
Símbolo de conformidad
Marque de conformité
Marchio di conformità

Not for reuse
Nicht zur Wiederverwendung
Ne reutilizar
Ne pas réutiliser
Non adatto per il riutilizzo

Note the instruction for use
Gebrauchsanweisung beachten
Consultar instrucciones para el uso
Observer le mode d'emploi
Osservare istruzioni per l'uso

Manufacturer
Hersteller
Fabricant
Fabricante
Fabbricante

**MODE D'EMPLOI
Bone Protect Fleece /
Bone Protect Cone**

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Bone Protect Fleece / Bone Protect Cone

Ref.	Size	Contents
Fleece - 12	2 cm x 2 cm	12 pieces
Cone - 12		12 pieces

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Bone Protect Fleece / Bone Protect Cone
Absorbable, Local Hemostat
Sterile

For single use only

Porcine collagen manufactured from porcine corium

100 g contain a minimum of 96,75 g collagen of porcine origine

Interaction with other agents

Antiseptics, which release chlorine (e.g. chloramine) as well tannic acid and caustic which change the protein must not be used together with collagen.

Ointments, powder as well silicon preparations can change intercostale space at collagenfoils, -sponge and -fleece, and should not be used together with collagen.

Indication areas

Bone Protect Fleece / Bone Protect Cone is indicated wherever capillary, venous, small arterial and diffuse seeping bleedings must be stopped, and where conventional means of hemorrhage control are either inadequate or technically difficult and time-consuming.

Indications**Maxillofacial Surgery**

- Socket preservation
- Hemostyptic treatment of extraction wounds of patients at risk of bleeding
- Sinus floor elevation
- Periimplant bone defects when placing immediate implants
- Filling the expansion gap in bone splitting
- Filling of bone defects (e.g. after cystectomy)

Dosage and Application

After opening of the package, Bone Protect Fleece / Bone Protect Cone is taken out with dry sterile instruments, placed on the cleaned wound and pressed on lightly with swabs. Bone Protect Fleece / Bone Protect Cone can be cut to the necessary size with sterile scissors. Immediately after contact with the wet wound surface Bone Protect Fleece / Bone

GB

Protect Cone sticks to the wound and forms a gel like bond with blood. With superficial wounds an additional layer of cotton gauze fixed with elastic bandages is recommended.

Contraindications

Bone Protect Fleece / Bone Protect Cone is contra-indicated in infected wounds, and must not be used concurrently with methyl methacrylate adhesive, it must not be used for intravascular application.

If you use an absorbable hemostatic agent on or near bony or neural spaces:

- use the minimum amount necessary to achieve hemostasis; and,
- remove as much of the agent as possible after hemostasis is achieved.

This will reduce the likelihood of neural and other soft tissue damage from swelling of the absorbable hemostatic agent, and/or migration and swelling of fragments of the agent.

Side effects

Allergic reactions may occur in isolated cases. Sometimes pain appears after application of a dry preparation to the wound surface. The increase of existing infection is very rare.

Note!**Bone Protect Fleece / Bone Protect Cone**

- is sterilized by γ -irradiation and must not be resterilized
- must be considered sterile and must be discarded if the package is damaged
- must not be used after the expiry date
- must not be stored above 24°C
- must be protected against moisture
- is a medical device and therefore must be kept from children!

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