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## CompactBone S

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## Instructions for use

## CompactBone S

Before using the product, please read these instructions for use carefully!

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

Note: CompactBone S is for professional use only. Responsibility for proper selection of patients, for adequate training, for experience in the choice of CompactBone S, all aspects of the surgery, and for the choice of post-operative procedures rests entirely with the surgeon.

## 1. Product description/ Material

CompactBone S is a mixture having a nominal composition of 60% hydroxyapatite and 40%  $\beta$ -tricalcium phosphate. These materials have been subjected to extensive and successful clinical studies for many years. CompactBone S is a safe and fully biocompatible material which is designed to act as an osteoconductive scaffold to support the ingrowth and fusion of adjacent viable bone when placed in an osseous environment. CompactBone S will later resorb and be replaced by natural bone.

CompactBone S can in many cases spare the patient the trauma of autograft harvesting and provides a sterile alternative to human or animal cadaver bone that is free of the risk of disease transmission.

## 2. Application areas

CompactBone S is a synthetic osteoconductive bone substitute in a fine granular form that is intended to be used as a scaffold to support the ingrowth of adjacent viable bone in defects that are not intrinsic to the stability of the bone structure. CompactBone S has been specially designed for filling of infrabony periodontal and root resection defects, filling of extraction sockets, elevation of the maxillary sinus floor, and alveolar ridge augmentation and reconstruction. CompactBone S may also be used in conjunction with products intended for guided bone regeneration.

The granules may be gently pressed into the defect and may also be mixed with sterile saline, blood, bone marrow aspirate or autogenous bone prior to implantation. CompactBone S will act as a temporary scaffold and is not intended to provide structural support during the healing process. The bone substitute material is biocompatible, radio-opaque and will eventually resorb in a controlled way.

## 3. Contraindications

CompactBone S is not designed or sold for any other use except as indicated.

Do not use CompactBone S in the presence of any contraindication.

CompactBone S is contraindicated where the device is intended to provide structural support and must not be used for the fixation of titanium fixtures.

Other contraindications include:

- Existing acute or chronic infections, particularly at the implantation site
- Severe vascular or neurological disease
- Uncontrolled diabetes
- Severe degenerative disease
- Hypercalcaemia, abnormal calcium metabolism
- Inflammatory bone disease
- Malignant tumours
- Severely impaired renal function

## 4. Precautions

It is important to maximise the contact between existing bone and the bone substitute material to ensure proper bone regeneration.

CompactBone S must not directly be used to achieve implant fixation. The implantation of titanium fixtures should not take place until about 6 months after the use of CompactBone S in an implant site.

## 5. Possible Complications

A successful result may not be achieved in every case.

A secondary operation to remove or replace the bone substitute material may be necessary due to surgical error, specific medical conditions or device failure. Possible adverse effects may include but are not limited to:

- Fracture or crushing of the bone substitute material with or without generation of particulate debris due to a load being applied
- Bone deformity or non-viable bone at the implantation site
- Post operative complications and other complications that are possible with surgery.
- Allergic reaction to the product, although no adverse reactions have been reported

## 6. Storage and Shelf Life

Store in a dry place at room temperature.

The expiration date is printed on the label. Do not use CompactBone S after the expiration date.

CompactBone S is environmentally friendly. No special disposal is required.

## 7. Application

Step 1: Expose the defect with mucoperiosteal flap, and carefully remove all granulation tissue

Open both the outer and inner pouches. Open the container.

Step 2: Implant the material. The bone substitute material may be first mixed with sterile saline, blood, bone marrow aspirate or autogenous bone as per normal surgical practice and may then be gently and carefully tamped into place. Excessive force will crush the trabecular structure. Overfilling of the defect should be avoided. The bone substitute material should only be placed in direct contact with well vascularised bone, and cortical bone should be mechanically perforated.

Step 3: Secure the surgical site after implanting to prevent implant migration. If possible, the mucoperiosteal flaps should be sutured to achieve primary closure. A surgical dressing may be placed over the wound for one to two weeks. Allow approximately six months for the bone to satisfactorily colonise the bone substitute material before placing titanium fixtures.

## 8. Warnings

The contents of the package are double sealed and STERILE. Do not use if opened, punctured or product appears damaged. Read the expiration date before use and do not use beyond the expiration date.

It is recommended not to open the inner pouch until a few minutes prior to implantation.

CompactBone S is for single use only. Do not attempt to re-sterilise or re-use due to the risk of cross infection.

The granules must be secured to prevent any potential migration of the bone substitute material and should only be used in procedures where the implant can be adequately contained.

CompactBone S is radiopaque and may hide areas under or above the bone substitute material on a radiograph.

## 9. Symbols

- Can be used until JJJJ-MM
- Lot number
- Sterilization by irradiation
- Do not reuse if packaging is damaged
- Single use
- Note: Observe accompanying documents

## 10. Other information

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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Dentegris Deutschland GmbH