BoneProtect® Guide

DE	GEBRAUCHSINFORMATION
EN	INSTRUCTIONS FOR USE
FR	INSTRUCTIONS D'UTILISATION
ES	INSTRUCCIONES DE USO
IT	ISTRUZIONI PER L'USO
PT	INFORMAÇÕES DE UTILIZAÇÃO
HU	HASZNÁLATI ÚTMUTATÓ
HR	INFORMACIJA O KORIŠTENJU
TR	KULLANIM BILGILERI
RU	ИНФОРМАЦИЯ ПО ПРИМЕНЕНИЮ
EL	ΠΛΗΡΟΦΟΡΙΕΣ ΧΡΗΣΗΣ

CE 0483

Fabricant / Manufacturer
MBP
Medical Biomaterial
Products GmbH
Lederstraße 7
D-19306 Neustadt-Glewe

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- Porcine collagen, produced from porcine corium
- For single use

The following sizes are available:

- 15 x 20 mm
- 20 x 30 mm
- 30 x 40 mm

Areas of application

Used solely or in combination with appropriate augmenting material (e.g. autogenic bone, allogeneic, xenogeneic or alloplastic bone replacement materials), BoneProtect® Guide is indicated for immediate or delayed controlled regeneration of tissue and bone in dental surgery.

- in the case of surgical bone and bone walls defects
- in connection with sinus floor augmentation and for supporting the Schneider's membrane
- in connection with alveolar ridge augmentation
- in connection with alveolar ridge reconstruction for prosthetic restoration
- in connection with treating fenestration defects
- in the case of periodontal bone defects (one to three-wall defects, functional defects category I, II)
- in the case of peri-implant dehiscence defects
- after root tip resection, removal of cysts, removal of retained teeth and resection of other bone defects
- within or on extraction sockets after tooth extraction
- in the case of immediate or delayed augmentation around implants within extraction sockets

Contraindications

- BoneProtect® Guide may not be used on patients suffering from
- acute infections of the oral cavity or acute or chronic infections of the implantation location
 - general illnesses which prohibit oral-maxillofacial surgical, implant, periodontal, endodontic or other oral surgery measures
 - known sensitivity to porcine collagen

Use during pregnancy and breastfeeding

There has neither been any testing on the use of BoneProtect® Guide during pregnancy and breastfeeding nor on its impact on human fertility. Therefore the advantages for the mother and the possible risks for the child must be assessed individually before using BoneProtect® Guide.

Use for children and older patients

There is no evidence that special precautions are necessary depending on the age of the patient to be treated.

Usage instructions

BoneProtect® Guide should only be applied by dentists who are trained in and acquainted with the techniques of controlled bone and tissue regeneration. BoneProtect® Guide is of natural origin. This is why a slight corrugation of the collagen structure and a variation of the membrane thickness may occur in the dry state. These occurrences do not have any effect on the quality or functionality of BoneProtect® Guide. BoneProtect® Guide can be adapted to the size and shape of the defect to be treated with a pair of scissors in a dried and rehydrated state, while maintaining sterility. The application of appropriate stencils can be useful for determining the required amount of BoneProtect® Guide. BoneProtect® Guide should overlap the walls of the defect by at least 2–3 mm. This enables total placement on the bone to prevent connective tissue from growing in on the sides. While using BoneProtect® Guide the general principles of sterile working and patient medication must be adhered to

- The required surgical treatment starts with presenting the defect
- The indicated bone defect is refilled with appropriate augmenting material (e.g. autogenic bone, allogeneic, xenogeneic or alloplastic bone replacement materials)
- The exterior packaging, which is sterile on the inside, must be removed e.g. by an assistant in its box in the unsterile operating area and opened while maintaining sterility. The interior packaging, which is sterile on the outside and inside, should be removed while maintaining sterility
- BoneProtect® Guide is placed on the defect and adapted to the spot with moderate pressure. The duration of pressure depends on the amount of bleeding. Adhesion to the bone surface is enabled

through swelling and gel formation of the collagen fibre in contact with blood. BoneProtect® Guide can be applied in a dried and wet state. If it is preferred to use the membrane in a wet state, the membrane should be rehydrated in a sterile saline solution before it is applied

- Complete soaking of BoneProtect® Guide with blood and exudate enables perfect adaptation and adhesion to the defect structure resp. the augmentation as well as development of a blood coagulum beneath the membrane
- After surgical treatment the patient should carry out oral hygiene measures according to the dentist's instructions

Special information on use in periodontology

Successful periodontal treatment is based on controlling the bacterial infection by means of debridement (removal of granulation tissue, subgingival curettage, scaling, root planning etc.), antibiotics treatment if applicable, and sufficient oral hygiene practiced by the patient according to the dentist's instructions. It is important to remove as little soft tissue as possible in order to enable ideal wound closure. Surgical treatment should be preceded by a hygiene phase with patient instruction and a renewed evaluation of the clinical situation by the dentist. In order to guarantee long-term treatment success, it should be followed by a post-operative maintenance phase with respective patient instruction by the dentist. To effectively prevent the development of a junctional epithelium, BoneProtect® Guide must be modelled closely to the tooth resp. the tooth neck and fixated with surgical suture, nails or pins, if necessary.

Dosage

The required amount of BoneProtect® Guide depends on the respective anatomic relations and the implant applied. BoneProtect® Guide is applied to the bone defect in the required amount and should overlap the walls of the defect by at least 2–3 mm. For this purpose BoneProtect® Guide can be cut to the right size with a pair of scissors. Suitable sterile stencils can help to determine the required size.

Side effects

- Allergic reactions to the collagen membrane cannot be excluded in rare cases
- There may be appearances of incompatibility with collagen in very rare individual cases
- There may be infectious tissue reactions due to a prolonged resorption time in very rare individual cases
- As with any foreign material, existing local infections can be intensified by the implantation of BoneProtect® Guide
- Possible general complications may result from the surgical treatment itself, such as recession of the gingiva, heavy bleeding of the gums, swellings of the soft tissue, sensitivity to temperature, desquamation of the gingiva epithelium in the area of the lobe

Interactions with other medication

- Antiseptics releasing chlorine (e.g. chloramine) as well as tannin and caustics, which modify proteins, may not be applied in combination with collagen
- Creams and powders as well as silicone compounds may modify the spaces between collagen foils, tiles or sponges, which is why they should not be used either

Warnings, precautions

BoneProtect® Guide is elastic and adheres to the bone. In order to create space and promote bone development, BoneProtect® Guide can be used with appropriate augmentation material (e.g. autogenic bone, allogeneic, xenogeneic or alloplastic bone replacement materials), if required.

The resorption time may be accelerated if BoneProtect® Guide is exposed during the healing phase.

BoneProtect® Guide is only indicated for the fields of application stated. The practitioner is responsible for informing the patients about possible contraindications, side effects and necessary precautions. In the case of post-operative complaints such as pain, infections or other unusual symptoms, the patient should seek immediate help from a dentist. Patients with severe general illnesses (e.g. badly adjusted diabetes mellitus, high blood pressure, peripheral artery occlusive disease, malignant or autoimmune disease) or patients requiring long-

term steroid treatment or anticoagulative therapy must be treated with special caution, just like in the case of any other surgical treatment.

Durability

The expiry date is printed on the exterior packaging and the sterile interior packaging. After expiry of the date stated on the packaging BoneProtect® Guide should no longer be used.

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